



The Fotoadapter allows you to set up a system for taking photos of microscopic specimens using Bresser MicroCam PC goggles and cameras from other manufacturers, provided that they are equipped with a C-mount (M25 thread, camera standard). Dedicated for use with scientific microscopes Bresser BioScience Trino and Bresser Science TRM-301. The photo adapter is placed on the photomicrograph stem. This method of assembly ensures simultaneous viewing and photographing of observed objects. Examples of cooperating cameras: Bresser Science, Bresser Microcam, Imaging Source, HDCE-50B 5MP. In the case of Bresser Microcam cameras, it is necessary to disassemble the camera element containing the 0.5x lens. (adapter + camera equipped with C-mount)