



The cross table is fixed by screwing it to the microscope table. It allows you to precisely position the desired fragment of the preparation in the field of view of the microscope, thanks to the possibility of moving it in two directions in the table plane. RANGE OF RUBBER: 30 mm in depth, 40 mm in width. Works with Bresser microscopes: BIOLUX NG, Erudit, Biolux AL, and VISIOMAR, as well as with many microscopes: Levenhuk 50L, 50L Plus, D50L, 2L, 2L Plus, selected Delta Optical, TPL, Ecotone etc. Warranty 2 years